

Chest Electrode (Child) with Logo

General

The chest electrode consists of a rubber suction bulb and a tip electrode connector which connects to an ECG electrode lead.

Use a tip type ECG electrode lead to connect the chest electrode to the electrocardiograph. Use an electrode lead with a tip size which fits the tip electrode connector.

Safety Information

WARNING A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

CAUTION A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

WARNING

- Do not use the chest electrode during MRI examination because it may cause skin burn on the electrode placement area.
- When performing defibrillation, discharge as far as possible from electrodes on the chest of the patient. If there is a possibility that the defibrillator paddle could touch electrodes, remove the electrodes from the patient. If the defibrillator paddle directly contacts an electrode, the discharged energy may cause skin burn to the patient.

CAUTION

- ECG measurement should be performed in as short a time as possible. If the chest electrode is used for a long time, it may cause redness on the skin.
- Do not immerse the chest electrode in water or saline solution for more than an hour. The electrode may get corroded or rusted and the ECG may be distorted.
- Store the chest electrode away from children.

CAUTION

Dispose of Nihon Kohden products according to your local laws and your facility's guidelines for waste disposal. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, dispose of it as medical waste according to your local laws and your facility's guidelines for medical waste. Otherwise, it may cause infection.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiuchi, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

EC REP European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599