

Interface

QI-670P

General

The QI-670P interface enables a BSM-6000 series bedside monitor to receive monitoring data from a TEC-5600 series or TEC-8300 series defibrillator via Bluetooth communication. The transmitted data can be saved in the QI-670P interface and reviewed on the bedside monitor.

Safety Information

- ⚠ WARNING** A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.
- ⚠ CAUTION** A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

⚠ WARNING

Never use the QI-670P interface in the presence of any flammable anesthetic gas or high concentration oxygen atmosphere. Failure to follow this warning may cause explosion or fire.

⚠ WARNING

Never use the QI-670P interface in a hyperbaric oxygen chamber. Failure to follow this warning may cause explosion or fire.

⚠ WARNING

Do not take the interface into the MRI test room. It is not designed to be used during MRI tests.

⚠ WARNING

Do not disassemble, modify or repair any parts that are not specified in this manual. Attach the cover to the interface after connecting it to the bedside monitor. Otherwise, the patient or operator may receive electrical shock. If you need to remove the cover for maintenance, disconnect the bedside monitor from the patient beforehand. Maintenance must be performed by a qualified service personnel.

⚠ CAUTION

Before connecting or disconnecting instruments, make sure that each instrument is turned off and the power cord is disconnected from the AC socket. Otherwise, the patient or operator may receive electrical shock or injury, data may be lost or the instrument may malfunction.

⚠ CAUTION

If fluids are accidentally spilled into the QI-670P interface, take the QI-670P interface out of service and contact your Nihon Kohden representative. The QI-670P interface must be disassembled, cleaned, dried and tested for safety and function.

⚠ CAUTION

Dispose of Nihon Kohden products according to your local laws and your facility's guidelines for waste disposal. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, dispose of it as medical waste according to your local laws and your facility's guidelines for medical waste. Otherwise, it may cause infection.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.



Manufacturer

NIHON KOHDEN CORPORATION
1-31-4 Nishiuchiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>



European Representative

NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599