

Battery Pack

SB-310V

General

The SB-310V battery pack is for an AED-3100 automated external defibrillator.

Safety Information

⚠ WARNING A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

⚠ CAUTION A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

⚠ WARNING

Never do the following. The following may cause leak, overheating, explosion and fire.

- Giving an impact to the battery pack by dropping or throwing.
- Charging, short-circuit, disassembly, deformation, overheating, dropping into fire, or immersing in water.

⚠ WARNING

The battery pack must be inserted in the AED when you install the AED. If there is no battery pack in the AED, self test is not performed and the AED condition is not maintained.

⚠ WARNING

If the battery pack is damaged and the substance inside the battery pack contacts the eyes or skin, wash immediately and thoroughly with water and see a physician. Never rub your eyes, because you may lose your eyesight.

⚠ WARNING

Do not expose the battery pack to direct sunlight or leave in a high temperature place such as in a car on a hot day or in front of a heater. The lifetime of the battery pack may be shortened, the performance of the battery pack may be degraded and the battery pack may leak.

⚠ WARNING

Do not use the battery pack if it is wet. Too much electricity might be discharged and the battery pack gets damaged.

⚠ WARNING

Do not use a broken or deformed battery pack. It may explode or cause fire.

⚠ CAUTION

Check the remaining battery power periodically.

⚠ CAUTION

Store the battery pack in the following conditions. Otherwise, the battery pack might deteriorate.

- Temperature: -20 to $+70^{\circ}\text{C}$ (-4 to 158°F)
- Humidity: 5 to 95% (noncondensing)

⚠ CAUTION

When disposing of the battery pack, follow your local law.

⚠ CAUTION

In the following cases, the life of the battery pack may become shorter than 4 years.

- Turning the AED on and off too frequently for purposes other than rescue or inspection.
- Turning the AED power on for a long time, such as during training or inspection.
- Discharging energy too many times for testing.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.

 Manufacturer
NIHON KOHDEN CORPORATION
 1-31-4 Nishiochiai, Shinjuku-ku,
 Tokyo 161-8560, Japan
 Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

EC REP European Representative
NIHON KOHDEN EUROPE GmbH
 Raiffenstrasse 10, D-61191 Rosbach, Germany
 Phone +49 6003-827-0 Fax +49 6003-827-599