INSTRUCTION FOR USE

Cleanac, Cleanac.3, Cleanac.710, Isotonac.4, Hemolynac.3N, Hemolynac.5, Hemolynac.310, Hemolynac.510.

IFU NUMBER: NKME/IFU-02

ISSUANCE DATE: 08th JUNE 2022

ENGLISH

INTENDED USE

For in vitro diagnostic use only. Reagents to be used on NIHON KOHDEN cell counters in combination as indicated in **Table 1**. For the intended use see **Table 2**. For professional use only.

METHODOLOGY AND CALCULATION OF THEANALYTICAL RESULT

Refer to the Operating Manual of the instrument.

STORAGE AND STABILITY AFTER FIRST OPENING

Store between 1 and 30 °C. Install each reagent after checking the integrity of the packaging and the validity of the device. Once opened, the reagent may be used until depleted and by the date indicated in **Table 3**. Use the undamaged original cap to reseal the packaging in the case of temporary cessation of use, except for Hemolynac.310, Hemolynac.510, Cleanac.710. Reseal without forcing.

COMPOSITION AND TARGET POPULATION

Active ingredients See **Table 4**. Target population see **Table 5**.

WARNINGS

Labeling information See **Table 6**. Any serious incidents that may occur involving the device must be reported both to the manufacturer and to the competent authority of the Member State in which the incident has occurred. **Information on disposal:** Product remnants, waste from its use and empty containers must be disposed of in accordance with applicable legislation. Waste from the instrument, such as potentially infected substances, must be handled with care and disposed of in accordance with the applicable legislation.

SAMPLE

The reagents were created for use with whole blood samples. Refer to the Operating Manual of the instrument for sample collection and storage procedures.

MARKER

Cleanac, Cleanac-3, Cleanac-710: N / A

Isotonac-4: Complete blood test in whole blood.

Hemolynac-3N, Hemolynac-310: Concentration of hemoglobin in whole blood.

Hemolynac-5, Hemolynac-510: Five differential white blood cells in whole blood.

PROCEDURE AND MEASUREMENT PRINCIPLES

Please refer to the Operating Manual of the instrument for information on how to install the reagents, how to preflush the lines, and how to use the instrument. Operate at the operating temperatures indicated in the Operating Manual of the instrument. Use each product in its original packaging (do not decant). Operate according to good laboratory practice (GLP).

Refer to **Table 7** for the measurement principles.

INTERFERING SUBSTANCES OR LIMITATIONS

Below the indicated concentrations, the following substances do not interfere:

Substance	Max. serum concentration			
	MEK-680D	MEK-910D	MK-310WD	MK-510WD
Bilirubin C	33.1 mg/dL	33.1 mg/dL	25.0 mg/dL	35.3 mg/dL
Bilirubin F	31.6 mg/dL	31.6 mg/dL	17.9 mg/dL	35.4 mg/dL
Hemolysis hemoglobin			0.18 g/dL	0.88 g/dL
Chyle	2384 FTU	2384 FTU	1031 FTU	2578 FTU
Total protein			1.02 g/dL	5.11 g/dL

INTERNAL QUALITY CONTROL

Background check

After the instrument lines have been preflushed, check the background (background count) as explained in the Operating Manual of the instrument, and make sure that the values obtained for each parameter are within the limits specified in the Operating Manual. If one or more of the parameters are not within these limits, proceed according to the instrument manufacturer's instructions. If the problem can be attributed to one of the reagents, replace the package in use with an intact one and repeat the check. Do not proceed with the analysis until the background values are acceptable. Calibration and calibration verification procedures. It is recommended to check the calibration of the instrument under the following conditions:

1. after the initial installation of the instrument
2. at least every six months (or, in the case of more restrictive indications, as recommended by the manufacturer of the instrument)
3. after a technical intervention deemed critical (replacement of components that may affect instrument performance, software upgrade, complete replacement of reagents)
4. if quality control data reveals a significant change in baseline accuracy (drift), or the data is outside the established acceptability limits.

Before performing calibration procedures, it is recommended to carry out preventive maintenance (baseline accuracy and precision check: ref. CLSI document H26-A2). Calibration of the instrument and calibration checking must be performed using appropriate material (fresh whole blood samples with values that have been assigned using manual reference methods or an appropriate calibrator), following the instructions provided in the Operating Manual of the instrument. If the instrument needs to be calibrated, it is recommended that the calibration check be performed immediately after calibration, following the instructions in the Operating Manual of the instrument (also refer to CLSI document H26-A2).

Daily quality control procedures

In order to ensure the accuracy and reliability of the results obtained, it is recommended to perform the quality control procedures following the instructions in the Operating Manual of the instrument, using appropriate control material (appropriate stabilized control blood), immediately before, during and immediately after each analysis sequence. It is recommended to use two different levels of stabilized control blood (refer to CLSI document H26-A2).

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IFU NUMBER: NKME/IFU-02
ISSUANCE DATE: 08th JUNE 2022

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INSTRUCTION FOR USE	
Cleanac, Cleanac.3, Cleanac.710, Isotonac.4, Hemolynac.3N, Hemolynac.5, Hemolynac.310, Hemolynac.510.	
IFU NUMBER: NKME/IFU-02	ISSUANCE DATE: 08 th JUNE 2022

Table 4 / الجدول 4 / Tableau 4 / Tablo 4 / Tabelle 4 / Tabella 4

Reagents / الكواشف / Réactifs / Reaktiver / Reagenzien / Reagenti		Active Ingredients / المكونات الفعالة / Ingrédients actifs / Aktif bileşenler / Wirkstoffe / Ingredienti attivi	
Name / الاسم / Nom / İsim / Name / Nome	REF	Name / الاسم / Nom / İsim / Name / Nome	Content / المحتوى / İçerik / Contenu / Inhalt / Contenuto
Isotonac-4	MEK-641D	Sodium chloride / كلوريد الصوديوم / Chlorure de sodium / Sodyum klorür / Natriumklorid / Sodio cloruro	4.0 g/L
		Anhydrous sodium sulfate / كبريتات الصوديوم اللامائية / Sulfate de sodium anhydre / Susuz sodyum sülfat / Wasserfreies Natriumsulfat / Sodio solfato anidro	9.8 g/L
Cleanac-3	MEK-620D	Sodium hypochlorite (% active Cl) / هيبوكلوريت الصوديوم (% كلور نشط) / Hypochlorite de sodium (% Cl actif) / Sodyum hipoklorit (% aktif Cl) / Natriumhypochlorit (% aktives Cl) / Ipoclorito di sodio (% Cl attivo)	1 - 1.5% p/v-w/v.
Cleanac	MEK-520D	Polyoxyethylene tridecyl ether / إيثر بولي أوكسي إيثيلين ثلاثي الدسييل / Éther tridécylique de polyoxyéthylène / Polioksietilen tridesil eter / Polyoxyethylentridecylether / Poliossietilene tridecil etere	0.5 g/L
Cleanac-710	MK-710WD	Sodium alkyl ether sulfate / كبريتات إيثر ألكيل الصوديوم / Sulfate d'éther alkylque de sodium / Sodyum alkil eter sülfat / Natriumalkylethersulfat / Sodio alchil etere solfato	0.7 g/L
Hemolynac-3N	MEK-680D	Quaternary ammonium salts / أملاح الأونيوم الرباعية / Sels d'ammonium quaternaire / Kuaterner amonyum tuzları / Quartäre Ammoniumsalze / Sali di ammonio quaternario	37,2 g/L
Hemolynac-310	MK-310W	Quaternary ammonium salts / أملاح الأونيوم الرباعية / Sels d'ammonium quaternaire / Kuaterner amonyum tuzları / Quartäre Ammoniumsalze / Sali di ammonio quaternario	27,0 g/L
Hemolynac-510 Hemolynac-5	MK-510WD MEK-910D	Sodium alkyl ether sulfate / كبريتات إيثر ألكيل الصوديوم / Sulfate d'éther alkylque de sodium / Sodyum alkil eter sülfat / Natriumalkylethersulfat / Sodio alchil etere solfato	1 g/L

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EN Detergent for Hematology analyzers mainly composed of surfactant. The biggest factors that pollute the flow path of a hematology analyzer are the proteins and lipids in the blood, which fall into the "hard" category of dirt. Therefore, this product is prepared for powerful removal of such strong blood stains. It contains a strong surfactant that solubilizes blood proteins and lipids, a metal chelating agent that serves as a cleaning aid, and a buffer. By flowing a constant amount of this detergent through the flow path of Hematology analyzer at regular intervals, dirt adhering to the flow path can be cleaned.

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EN Detergent for Hematology analyzers mainly composed of hypochlorous acid. The biggest factors that pollute the flow path of a hematology analyzer are the proteins and lipids in the blood, which fall into the "hard" category of dirt. Therefore, this product is prepared for powerful removal of such strong blood stains. It contains a hypochlorous acid that chemically dissolves blood proteins and lipids. By flowing a constant amount of this detergent through the flow path of Hematology analyzer at regular intervals, dirt adhering to the flow path can be cleaned.

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Isotonac.4

MEK-641D

EN Diluent for measuring total blood cell count (CBC). It is a conductive isotonic solution with a pH close to that of blood, so the morphology of diluted blood cells is maintained during measurement.

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EN Lysing reagent: Hemolynac.3N is a hemolytic reagent for hemoglobin Measurement. Hemolynac.3N lyses the red blood cells to elute hemoglobin in the diluted blood sample so that the hemoglobin concentration can be measured. The eluted hemoglobin reacts with the quaternary ammonium salt in the reagent and changes to a hemoglobin compound. The hemoglobin compound is measured by absorbance (540-555 nm) by the Nihon Kohden family of hematology analyzers. The absorbance of the hemoglobin compound is proportional to the hemoglobin concentration. This hemolytic reagent does not contain cyanide.

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IFU NUMBER: NKME/IFU-02

ISSUANCE DATE: 08th JUNE 2022

EN Lysing reagent: Hemolynac-310 is a hemolytic reagent for hemoglobin Measurement. Hemolynac-310 lyses the red blood cells to elute hemoglobin in the diluted blood sample so that the hemoglobin concentration can be measured. The eluted hemoglobin reacts with the quaternary ammonium salt in the reagent and changes to a hemoglobin compound. The hemoglobin compound is measured by absorbance (520 nm) by the Nihon Kohden family of hematology analyzers. The absorbance of the hemoglobin compound is proportional to the hemoglobin concentration. This hemolytic reagent does not contain cyanide.

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EN Lysing reagent: it is a hemolytic reagent for the five differential white blood cells measurement. It lyses the red blood cells to measure the five differential white blood cells in the diluted blood sample so that Lymphocytes, Monocytes, Neutrophils, Eosinophils and Basophils can be measured. The five differential white blood cells are measured by the flow cytometry method using a semiconductor laser of hematology analyzers manufactured by Nihon Kohden.

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