HEMATOLOGY CONTROL+CRP

MK-3CL, MK-3CN, MK-3CH

English

General

Intended Purpose

For in vitro diagnostic use.

HEMATOLOGY CONTROL+CRP is for in vitro diagnostic use as a value assigned whole blood control designed to monitor the values of Nihon Kohden hematology and clinical chemistry analyzers. HEMATOLOGY CONTROL+CRP is not used for diagnosis or assistance in diagnosis, prognosis, prediction of disease or patient's physiological condition. Measuring this HEMATOLOGY CONTROL+CRP ensures that the device is functioning according to its intended use.

Control parameters are below:

- 1) WBC: White Blood Cell Count
- 2) RBC: Red Blood Cell Count
- 3) HGB: Hemoglobin Concentration
- 4) HCT: Hematocrit Percent
- 5) MCV: Mean Corpuscular Volume
- 6) MCH: Mean Corpuscular Hemoglobin
- 7) MCHC: Mean Corpuscular Hemoglobin Concentration
- 8) PLT: Platelet Count
- 9) MPV: Mean Platelet Volume
- 10) RDW-CV: Red Blood Cell Distribution Width in Coefficient of Variation
- 11) LY%: Lymphocyte Percent
- 12) MO%: Monocyte Percent
- 13) GR%: Granulocyte Percent
- 14) LY: Lymphocyte
- 15) MO: Monocyte
- 16) GR: Granulocyte
- 17) CRP: C-reactive protein

Before and during use, read the analyzer operator's manual and the assay sheet (which is provided with HEMATOLOGY CONTROL+CRP) together with this manual.

- NOTE**
- Overall judgement must be performed by the physician, referring to the analysis result, clinical findings and other examination results.
 - Use HEMATOLOGY CONTROL+CRP with the specified analyzers only.

- Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative.

Devices Intended for Use in Combination with

MEK-1303 automated hematology and clinical chemistry analyzer

Materials Provided and Materials Required

Materials Provided

- MK-3CL HEMATOLOGY CONTROL+CRP
- MK-3CN HEMATOLOGY CONTROL+CRP
- MK-3CH HEMATOLOGY CONTROL+CRP

Materials Required (Not Provided with the Control)






- Protective gloves
- Protective glasses
- Waste container for biohazardous materials

Intended Users

For laboratory professional use only, in laboratories with suitable equipment for hematological testing. Qualified personnel, e.g. laboratory technicians trained in hematology analysis techniques, will be able to use according to the operator's manual.

Symbols

The following symbols are used with the control. The descriptions of each symbol are given in the table below.

Symbol	Description	Symbol	Description
	Caution		Temperature limits
	Biological risks		This way up
	In vitro diagnostic medical device		Open vial stability
	European representative		Unique Device Identifier
	Manufacturer		Control/Concentration level
	Lot number		
	Catalogue number		
	Use by		
	Operator's manual; operating instructions		

0614-908003

Symbol	Description
	<p>The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.</p>

Safety Information

⚠ WARNING A warning alerts the user to the possible injury or death associated with the use or misuse of the instrument.

⚠ CAUTION A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in this operator's manual.

⚠ WARNING

- **POTENTIALLY BIOHAZARDOUS MATERIAL.** For in vitro diagnostic use. Each human donor/unit used in the preparation of HEMATOLOGY CONTROL+CRP has been tested by a FDA licensed method/test and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of vials follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures.
- When handling HEMATOLOGY CONTROL+CRP, wear rubber gloves to protect yourself from infection.
- Do not swallow HEMATOLOGY CONTROL+CRP. It may be infectious.

- NOTE**
- Do not use the product if it is past the expiration date on the package or if the product is stored under unspecified conditions.
 - Store the control in a pharmaceutical refrigerator.
 - Store the control at temperatures between 2 and 8°C (36 and 46°F).
 - Do not freeze the control.
 - If the cap of the control becomes dirty, wipe with a non-woven towel such as Kimwipes.
 - Do not use the control if deterioration is suspected. After mixing, the control should be similar in appearance to fresh whole blood. In unmixed tubes/vials, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration.

Using the Control

NOTE: Refer to the analyzer operator's manual for details on measurement.

Measurement of a control is recommended for statistical quality control of the measurement accuracy of the analyzer. Also, it is recommended to use the average and upper and lower limits set by each laboratory to measure the control.

Frequency of Performing Quality Control and Restrictions on Performing Quality Control

Determine according to the procedures established by each laboratory.

Control Used

MK-3CL, MK-3CN, MK-3CH HEMATOLOGY CONTROL+CRP

Measurement Principle

Refer to the analyzer operator's manual.

Procedure

1. Take the control out of the pharmaceutical refrigerator and allow it to warm to room temperature (15 to 30°C, 59 to 86°F) for 15 minutes before mixing.

2. Mix the control.


NOTE

- Do not use a mechanical mixer.
- Do not mix the control roughly. If it bubbles, hemolysis is caused.

- 1) Hold the tube containing the control horizontally between the palms of the hands and roll it back and forth for 20 to 30 seconds. Occasionally invert the tube so that the control is mixed thoroughly.

Roll the tube back and forth

Invert the tube

- 2) Continue to mix in the manner described in step 1) until the red cells are completely suspended. Controls stored for a long time may require extra mixing.
3. Gently invert the tube until there is no clot of blood on the bottom before measurement.
4. Measure the control as instructed in the analyzer operator's manual.

5. If the tube was opened for the measurement, clean residual material from the cap and tube rim with a lint-free tissue after the measurement, then close the cap tightly.
6. Return the tube to the refrigerator within 30 minutes.

Technical Information

Metrological Traceability of Values Assigned to the Control

For each parameter, the internationally recognized reference method is:

- 1) WBC and RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3 Vol.20 No.28
- 3) HCT: CLSI H07-A3 Vol.20 No.18
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001
- 5) LY, MO and GR: CLSI H26-A2:Vol.29 No.40 and CLSI H20-A2:Vol.27 No.4
- 6) CRP: ERM-DA474/IFCC

Expected Values

- When entering the assay value into the analyzer, read all the three barcodes on the assay sheet.
- If the barcode is difficult to scan, place the assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2 cm.
- The CRP value on the assay sheet is a reference range. Accuracy cannot be monitored by the CRP value of HEMATOLOGY CONTROL+CRP. Use the CRP value for precision (reproducibility) monitoring.
- Verify that the lot number on the tube matches the lot number on the assay sheet.
- Refer to the assay values for the specific instrument model on the provided assay sheet.
- Assay values are determined on well-maintained, properly calibrated instruments using the conformed reagents. Reagent differences, maintenance, operation technique, and calibration may contribute to inter-laboratory variation.

Composition

The control contains human erythrocytes, simulated leukocytes, mammalian platelets and CRP suspended in a plasma-like fluid with preservatives.

Sterilization Method

The control is not intended to be sterilized or kept in a sterile environment.

Interfering Substances or Limitations

The performance of the control is assured only if it is properly stored and used as described in this operator's manual.

Incomplete mixing of the tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

Analytical Performance Characteristics

- Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.
- Assay values on a new lot of control should be confirmed before the new lot is put into routine use. Test the new lot when the instrument in good working order and quality control results on the old lot are acceptable. The laboratory's recovered mean should be within the assay range.
- For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range. The user may establish assay values not listed on the assay sheet, if the control is suitable for the method.

Mathematical Approach upon which the Calculation of the Analytical Result is Made

The following parameters are calculated as follows ¹:

$$\text{MCV} = \text{HCT} / \text{RBC}$$

$$\text{MCH} = \text{HGB} / \text{RBC}$$

$$\text{MCHC} = \text{HGB} / \text{HCT}$$

$$\text{LY}\% = \text{LY} / \text{WBC}$$

$$\text{MO}\% = \text{MO} / \text{WBC}$$

$$\text{GR}\% = \text{GR} / \text{WBC}$$

¹ Without unit correction

Environmental Conditions

Storage and Transport Environment

Temperature: 2 to 8°C (36 to 46°F)

Expiration Date

The expiration date is shown on the package.

Lifetime: 3 months

Shelf Life Date After Opening

Use the control within 14 days of opening.

Package and Catalog Number

Model	Qty	Catalog Number
MK-3CN	MK-3CN (2.5 to 2.9 mL) × 3	3CN03
	MK-3CN (2.5 to 2.9 mL) × 6	3CN06
	MK-3CN (2.5 to 2.9 mL) × 12	3CN12
MK-3CL, MK-3CH	MK-3CL (2.5 to 2.9 mL) × 3, MK-3CH (2.5 to 2.9 mL) × 3	3CLH3
	MK-3CL (2.5 to 2.9 mL) × 6, MK-3CH (2.5 to 2.9 mL) × 6	3CLH6

Disposal

WARNING

Dispose of HEMATOLOGY CONTROL+CRP according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for disposing of infectious medical waste. Otherwise, it may affect the environment. If there is a possibility that HEMATOLOGY CONTROL+CRP may have been contaminated with infection, it may cause infection.

When disposing of the control, such as when the expiration date is past, follow the instructions on the SDS of the control.

Revision History

Edition	Date	Details	Code Number
1st Edition	20 Nov 2017	Initial issue	0604-912978
5th Edition	08 Nov 2022	IVDR compliance	0614-908003

- NOTE**
- The code number of this manual was changed from 0604-912978C to 0614-908003 when the manual was updated from 4th Edition to 5th Edition.
 - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Note for users in the territory of the EEA and Switzerland:
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599
UK Responsible Person
NIHON KOHDEN UK LTD.
Unit 3, Heyworth Business Park,
Old Portsmouth Road, Peasmarsh,
Guildford, Surrey, GU3 1AF, UK
Phone +44 14-8333-1328

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- 1) Leucocytes et érythrocytes :
ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB : CLSI H15-A3 Vol.20 N° 28
- 3) HCT : CLSI H07-A3 Vol.20 N° 18
- 4) PLT : International Council for Standardization in Hematology Expert Panel on Cytometry et International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am. Journal of Clinical Pathology 115:460-464 2001
- 5) LY, MO et GR : CLSI H26-A2:Vol.29 No.40 et CLSI H20-A2:Vol.27 No.4
- 6) CRP : ERM-DA474/IFCC

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- 1) WBC und RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3 Vol. 20 Nr. 28
- 3) HCT: CLSI H07-A3 Vol. 20 Nr. 18
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001
- 5) LY, MO und GR: CLSI H26-A2:Vol. 29 Nr. 40 und CLSI H20-A2:Vol. 27 Nr. 4
- 6) CRP: ERM-DA474/IFCC

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- 1) WBC e RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol. 1988; 10:203-212
- 2) HGB: CLSI H15-A3 Vol.20 No.28
- 3) HCT: CLSI H07-A3 Vol.20 No.18
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001
- 5) LY, MO e GR: CLSI H26-A2:Vol.29 No.40 and CLSI H20-A2:Vol.27 No.4
- 6) CRP: ERM-DA474/IFCC

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