


**MEK-5DN** HEMATOLOGY CONTROL  
for WBC 5 part differential 3.0mL

**ASSAY SHEET**

LOT	MEK5D23YN
	2024-02-05
OPEN	14DAYS

MEK-9100,9200  
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 <sup>3</sup> /μL	±0.8	8.0	7.9	7.8	7.92
NE%	%	±15.0	46.0	46.0	46.0	45.56
LY%	%	±15.0	34.6	34.6	34.6	37.44
MO%	%	±5.0	11.6	11.6	11.6	8.94
EO%	%	±5.0	5.9	5.9	5.9	7.02
BA%	%	±4.0	1.9	1.9	1.9	1.04
NE	10 <sup>3</sup> /μL	±1.5	3.6	3.6	3.6	3.60
LY	10 <sup>3</sup> /μL	±1.5	2.8	2.7	2.7	2.97
MO	10 <sup>3</sup> /μL	±0.5	0.9	0.9	0.9	0.71
EO	10 <sup>3</sup> /μL	±0.5	0.5	0.5	0.5	0.56
BA	10 <sup>3</sup> /μL	±0.4	0.2	0.2	0.1	0.08
RBC	10 <sup>6</sup> /μL	±0.20	4.64	4.57	4.62	4.74
HGB	g/dL	±0.5	13.7	13.6	13.5	13.79
HCT	%	±3.0	42.7	40.9	40.2	47.5
MCV	fL	±6.0	92.0	89.5	87.0	100.3
MCH	pg	±2.8	29.5	29.8	29.2	29.1
MCHC	g/dL	±3.0	32.1	33.3	33.6	29.0
RDW-CV *	%	±4.0	14.4	15.5	14.6	15.3
PLT	10 <sup>3</sup> /μL	±50	265	282	259	272.7
MPV	fL	±2.0	3.3	5.1	3.5	6.8

Manual Registration	
Lot	23Y
Level	Normal
Control	MEK5D

\* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901691M