


MK-RE1 HEMATOLOGY CONTROL
 for reticulocyte 3.0mL

ASSAY SHEET

LOT	MKRE243-1
	2024-05-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	1.46	±1.54
RET	10 ⁶ /μL	0.0527	±0.0557

2D Barcode Entry



• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
 It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration						
Lot	243					
Level	1					
Expiration	5	M	5	D	2024	Y
Control	MKRE					



NIHON KOHDEN CORPORATION
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN