


**MEK-5DH** HEMATOLOGY CONTROL  
for WBC 5 part differential 3.0mL

**ASSAY SHEET**

LOT	MEK5D245H
	2024-08-05
OPEN	14DAYS

MEK-9100,9200  
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 <sup>3</sup> /μL	±2.2	25.3	25.2	24.8	26.13
NE%	%	±10.0	68.1	68.1	68.1	68.05
LY%	%	±10.0	14.7	14.7	14.7	16.67
MO%	%	±5.0	11.1	11.1	11.1	9.16
EO%	%	±5.0	4.7	4.7	4.7	4.81
BA%	%	±4.0	1.4	1.4	1.4	1.31
NE	10 <sup>3</sup> /μL	±2.8	17.2	17.1	16.9	17.78
LY	10 <sup>3</sup> /μL	±2.8	3.7	3.7	3.6	4.36
MO	10 <sup>3</sup> /μL	±1.0	2.8	2.8	2.8	2.39
EO	10 <sup>3</sup> /μL	±1.0	1.2	1.2	1.2	1.26
BA	10 <sup>3</sup> /μL	±0.8	0.4	0.4	0.3	0.34
RBC	10 <sup>6</sup> /μL	±0.25	5.58	5.54	5.53	5.80
HGB	g/dL	±0.7	17.9	17.6	17.7	17.91
HCT	%	±4.0	55.1	53.5	52.0	63.2
MCV	fL	±8.0	98.8	96.6	94.1	109.0
MCH	pg	±2.8	32.1	31.8	32.0	30.9
MCHC	g/dL	±3.0	32.5	32.9	34.0	28.3
RDW-CV *	%	±4.0	13.9	14.8	14.1	15.5
PLT	10 <sup>3</sup> /μL	±70	556	547	554	559.3
MPV	fL	±4.0	6.3	7.3	6.4	9.1

Manual Registration	
Lot	245
Level	High
Control	MEK5D

\* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-9016920