


<b>MEK-3DN</b>	HEMATOLOGY CONTROL	2.0mL
----------------	--------------------	-------

**ASSAY SHEET**

LOT	B248N
	2024-11-05
OPEN	14DAYS

<b>MEK-1301,1302,1305 Barcode Entry</b>
---

The barcodes may be read in any order.  
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
<b>WBC</b>	10 <sup>3</sup> /μL	±0.8	8.0	8.0	7.90
<b>RBC</b>	10 <sup>6</sup> /μL	±0.20	4.83	4.83	4.83
<b>HGB</b>	g/dL	±0.5	14.0	14.2	14.27
<b>HCT</b>	%	±3.0	41.6	41.6	47.6
<b>MCV</b>	fL	±7.0	86.1	86.1	98.5
<b>MCH</b>	pg	±2.4	29.0	29.4	29.5
<b>MCHC</b>	g/dL	±3.0	33.7	34.1	30.0
<b>PLT</b>	10 <sup>3</sup> /μL	±50	249	249	268.7
<b>MPV</b>	fL	±2.0	5.9	8.2	8.5
<b>RDW-CV *</b>	%	±3.0	16.6	16.6	16.4
<b>LY%</b>	%	±10.0	32.9	32.9	30.11
<b>MO%</b>	%	±5.0	3.7	3.7	3.31
<b>GR%</b>	%	±10.0	63.4	63.4	66.58
<b>LY</b>	10 <sup>3</sup> /μL	±1.0	2.6	2.6	2.38
<b>MO</b>	10 <sup>3</sup> /μL	±0.4	0.3	0.3	0.26
<b>GR</b>	10 <sup>3</sup> /μL	±1.0	5.1	5.1	5.26

<b>MEK-1305</b>
QC ESR
8±5 mm

Manual Registration	
Lot	248
Level	Normal
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



**NIHON KOHDEN CORPORATION**  
 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN



0694-901685Q