

Fighting Disease with Electronics



HEMATOLOGY REAGENT CATALOGUE

*Accurate diagnosis, trusted health-**Blood don't lie***




HEMATOLOGY REAGENTS

NIHON KOHDEN CORPORATE VALUES


Integrity

We act with high ethical standards to save the lives of many.



Customer Centric

We take every action to improve the experience of patients and customers.




Humbleness

We approach issues with an open mind and seek out objective facts.



Goal Oriented

We focus on reaching our destination and explore possibilities to get there.



Diversity

We respect and leverage individual strengths to maximize our potential.



Initiative

We communicate and contribute proactively to meet challenges.



Creativity

We imagine new ideas and solutions to create unprecedented value.




ABOUT THE REAGENT FACTORY

Hematology Reagents are hydrous solutions for diluting and/or treating blood samples. Nihon Kohden Celltac reagents are manufactured in our reagent Plant in Dubai, DAFZA. This plant was established in the year 2020.


The reagents manufactured by Nihon Kohden are made of highly purified freshwater of high-grade chemicals. Finally, Nihon Kohden's reagents meet the Ministry of Health and ISO requirements hereby enabling us to achieve one of the highest customer satisfaction in the Middle East and Africa.

KEY FEATURES AND BENEFITS




MANUFACTURED IN A HIGH-QUALITY FACTORY

That means ISO-9001 and ISO-13485 standards.




PERFORMANCE TESTED

Whilst strict specifications and large lot sizes ensure a minimal lot-to-lot variation.




REGISTERED AND CE CERTIFIED

Guaranteed that products passed the governing registration agencies.




LONG SHELF LIFE

All reagents have an un-open shelf life of up to 18 months and have extra strong storage boxes to prevent collapse of heavy diluent.



ENVIRONMENTAL PROTECTION

All reagents are cyanide-free, protecting animal, human and environmental health.



HIGH TECHNOLOGY

Our technology guarantees excellent product stability.

OUR PRODUCTS

CLEANAC 810 (MK-810W)



Description: Detergent

Volume: 15 ml x 3 bottles

Applicable Analyzer: MEK-9100, MEK-9200

Commodity: Dangerous Cargo UN1791, Class 8

Mode of Shipment: AIR / SEA / ROAD

(temperature controlled 2-8 degrees)



HEMOLYNAC 510 (MK-510WD)




Description: Lysing Reagent

Volume: 250ML

Applicable Analyzer: MEK-9100, MEK-9200

Commodity: General Cargo

Mode of Shipment: AIR / SEA / ROAD



HEMOLYNAC 310 (MK-310WD)



Description: Lysing Reagent

Volume: 250ML

Applicable Analyzer: MEK-9100, MEK-9200, MEK-1301, 1302, 1303, 1305

Commodity: General Cargo

Mode of Shipment: AIR / SEA / ROAD



RETICULONAC (MK-110W)



Description: Staining reagent (Dye acridine orange)


Volume: 250ML

Applicable Analyzer: MEK-9200

Commodity: General Cargo

Mode of Shipment: AIR / SEA / ROAD

(temperature controlled 2-25 degrees)



OUR PRODUCTS

CLEANAC 3 (MEK-620D)



Description: Detergent



Volume: 1L




Applicable Analyzer: MEK-6500, 6510, 6550, MEK-7300,
MEK-1301, 1302, 1303, 1305



Commodity: Dangerous Cargo UN1791, Class 8



Mode of Shipment: AIR / SEA / ROAD



HEMOLYNAC 5 (MEK-910D)



Description: Lysing Reagent for 5-Part Diff



Volume: 1L



Applicable Analyzer: MEK-7300



Commodity: General Cargo



Mode of Shipment: AIR / SEA / ROAD



HEMOLYNAC 3N (MEK-680D)



Description: Cyanide-Free Lysing Reagent for 3-Part Diff



Volume: 1L



Applicable Analyzer: MEK-6500, 6510, 6550, MEK-7300



Commodity: General Cargo



Mode of Shipment: AIR / SEA / ROAD



OUR PRODUCTS

CLEANAC 710 (MK-710WD)



Description: Detergent



Volume: 3L



Applicable Analyzer: MEK-9100, MEK-9200,
MEK-1301, 1302, 1303, 1305



Commodity: General Cargo



Mode of Shipment: AIR / SEA / ROAD



CLEANAC (MEK-520D)



Description: Detergent



Volume: 5L



Applicable Analyzer: MEK-6500, 6510, 6550, MEK-7300



Commodity: General Cargo



Mode of Shipment: AIR / SEA / ROAD



ISOTONAC 4 (MEK-641D)



Description: Diluent



Volume: 20L



Applicable Analyzer: MEK-1301, 1302, 1303, 1305, MEK-6500, 6510,
MEK-7300, MEK-9100, MEK-9200



Commodity: General Cargo



Mode of Shipment: AIR / SEA / ROAD



Fighting Disease with Electronics



CONTACT US

Nihon Kohden Middle East

 JAFZA ONE TOWER A Jebel Ali Freezone - Dubai, Office 1912

 <https://ae.nihonkohden.com/>

 +971 4 884 0080

 ivdsales@nkme.ae

LinkedIn



 SCAN ME

