

# ECG Electrodes for the NICU

## Designed for High-Humidity Incubators

**G300D (N-01IS3)**



- Durable and adhesive gel hardly swells, preventing electrodes from peeling even under humid environment
- Smallest electrodes among Vitrode series
- Soft radiolucent lead wires are bundled to prevent the leads being tangled around the baby

Target: low-birth weight neonate, infant

Length: 0.6 m

Quantity per box: 3 colors/pack \* 10 packs

## Skin-Friendly with Gentle Adhesive

**G300A (N-03IS3)**



- Skin-friendly adhesive gel; gentle for less skin irritation
- Electrodes for the baby whose dermis is immature
- Soft radiolucent lead wires are bundled to prevent the leads being tangled around the baby

Target: low-birth weight neonate, infant

Length: 0.6 m

Quantity per box: 3 colors/pack \* 10 packs

## Available from Neonate to Infant

**G271A (V-120S3)**



- Bitter chemical compound is used on the electrodes to prevent accidental swallowing
- Electrodes are specifically designed for the more active baby
- Soft radiolucent lead wires are bundled to prevent the leads being tangled around the baby

Target: neonate, infant




Length: 0.6 m

Quantity per box: 3 colors/pack \* 40 packs

# ECG Electrodes for the NICU



## Comparison Chart

Order Code (Product Code)	Size (mm)	Target	Adhesion		Skin Friendliness	Adhesive Strength to Body Movement
			Humidified Environment*	Non-Humidified Environment		
G300D (N-01IS3) 	15*15	Low birth weight neonate, infant	★ ★ ★ ★ ★	★ ★ ★	★ ★ ★ ★	★ ★
G300A (N-03IS3) 	14*25	Low birth weight neonate, infant	★	★ ★	★ ★ ★ ★ ★	★
G271A (V-120S3) 	20*20	Neonate, infant	★ ★	★ ★ ★ ★	★ ★	★ ★ ★ ★

\*Humidified Environment means inside incubator; indicates durability of electrodes in incubator.

• Star indicates the strength in each category. The more colored stars (★) the better the electrodes address the category.