

Bedside Monitor

SVM-7600 series



For advanced monitoring

Fighting Disease with Electronics

 **NIHON KOHDEN**

For advanced monitoring

Intuitive touch screen operation

The touch screen provides easy and intuitive operation. Just click the data area for the patient or parameter to see more information or change a setting. Easy operation frees you to focus on your patient.


120 hours of synchronized review data

Up to 120 hours of review data is saved.

Time is synchronized across all review screens so you can easily review an event different points of view. Comprehensive review contributes to decision making for patients.



Trend table




Full disclosure



Alarm history

Interbed gives quick access to check any patient in the network

From any monitor, you can check the vital information and alarm status of any other monitor in the network. With this interbed function, you can immediately check the status of any patient. This leads to better care for all patients.



Multiple Beds Window



Individual Bed Window

SVM-7603
10.4-inch LCD



SVM-7623
12.1-inch LCD




User-friendly operation

An onscreen guide with helpful illustrations makes it easy to operate the monitor and perform basic tasks even for first time users.

The guide also includes an explanation of abbreviations for parameters, arrhythmias and icons which are found in all Nihon Kohden monitors.

Smooth operation enables you to focus on better patient care.

Onscreen guide



Abbreviation list

PARAMETERS	ITEMS	EXPLANATION
CO ₂ -RR		Respiration rate (CO ₂)
DIA		Diastolic pressure
HR		Heart rate
IMP-RR		Respiration rate (impedance)
MAP		Mean arterial pressure
MEAN		Mean pressure
NEP		Non-invasive Blood pressure

Parameters

PARAMETERS	ITEMS	EXPLANATION
ASYSTOLE		Asystole
BRIDGEMAN		Ventricular bigeminy
BRADYCARDIA		Bradycardia
COUPLET		VPC couplet
EARLY VPC		Early VPC, including R on T type
FREQ VPC		Frequent VPC (reaching or exceeding the preset limit per minute)
TACHYCARDIA		Tachycardia

Arrhythmias


PARAMETERS	ITEMS	EXPLANATION
	<input type="checkbox"/>	All alarm off
	<input type="checkbox"/>	Alarm suspended
	<input type="checkbox"/>	Alarm silence
	<input type="checkbox"/>	Alarm silence with rearming time
	<input type="checkbox"/>	Alarm reset
	<input type="checkbox"/>	Battery status
	<input type="checkbox"/>	Checking data

ITEMS

Reliability, quality and performance

To assure high quality and performance, SVM-7600 bedside monitors are designed and manufactured under the same rigorous standards as all Nihon Kohden products.

Nihon Kohden's assured quality over the life of the product gives you excellent total cost of ownership and best return on investment.



This brochure may be revised or replaced by Nihon Kohden at any time without notice.



NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100
www.nihonkohden.com